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MANUEL D'INSTRUCTION

SOMMAIRE

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INSRTUCTION MANUAL

CONTENTS

FOREWORD	34
1. Safety Guide	34
1.1 Intended Use Of 4 Light	34
1.2 Checking The Mask	34
2. Technical Data	34
3. Device Parts	35
4. The LED Light Therapy Mask Description	35
4.1 The LED Light Therapy Mask Consists Of	35
4.2 Features	35
4.3 Remote Wireless Controller Marking	35
5. Contraindications, Warning And Precautions	36
5.1 Contraindications	36
5.2 Warnings	37
5.3 Precautions	38
6. Operating Mode Information	40
7. Sensitization Test	41
8. Using LED Light Therapy Mask	41
8.1 Treatment Recommendations	41
8.2 Carrying Outa Treatment	42
8.3 Tips	43

9. Using The Wireless Controller	43
9.1 Modes	43
9.2 Light Level Selection	43
9.3 Voice Announcement	43
10. Charging	44
11. Risks And Benefits	48
12. Safety Signs And Meaning	48
13. Care And Maintenance	49
13.1 Cleaning Of LED Light Therapy Mask	49
13.2 Disposal	49
14. Troubleshooting	50
15. Environmental Conditions	51
16. EMC Manufacturer's Declaration And Guidance	52
17. Service Life	58
18. Guarantee Information	58
19. Frequently Asked Questions	59

FOREWORD

This LED light therapy mask adopts the latest LED light therapy technology, combining 4 kinds of light, one machine can solve various skin problems.

LED light mask is a home use wearable LED device which can only be used in these scenarios and can not be used for other purposes, otherwise it may cause damage to the product or the human body. Manufacturers and distributors are not liable for damage caused by improper use.

The product is not used to make measurements of any sorter to draw any conclusions regarding the indication to treat. The equipment does not require checks on the light output as the LED do not dim with age to any practical extent.

Version identifier: 240924.

1.Safety Guide

The LED light therapy mask is an electrical device, which under certain circumstances could present an electrical shock hazard to the user. Please read the accompanying user manual thoroughly and follow directions stated in the manual to assure maximum safety during operation. The intended operator or user is in good health for use.

1.1 Intended Use Of 4 Lights

- (1) Red light: Assist in reducing wrinkles.
- (2) Deep red light: Assist in reducing deep wrinkles.
- (3) Blue light:Help decreasing mild to moderate inflammatory acne.
- (4) Infrared light: provide topical heating for the purpose of elevating tissue temperature; arthritis and muscle spasm; relieving stiffness; promoting the relaxation of muscle tissue; and to temporarily increase local blood circulation.


1.2 Checking The Mask

After receiving the product, if you use scissors or a knife to open the package, be careful not to scratch the product. In addition, it is necessary to check whether the product accessories are complete, whether the product can be used normally. If there is any omission or the product is not working properly, you should contact customer service to replace it.

2. Technical Data

Product Name: 4D GLOW LED LIGHT THERAPY MASK	
Model: FM19	Rated Voltage: 3.7V
Battery Capacity: 850mAh	Battery Cycle Life: ≥300times
Input Voltage : 100 V - 240 V (50/60 Hz)* See FAQ	
Charging Current : 5 V / 1 A (Recommended Amperage)	
Quantity of color LED lights: 60 pieces	
Dimension(L×W×H): 578×282×133 mm	
IP classification: IP22	Gross Weight: 286 g
The classification of LS EQUIPMENT: Exempt Group	

3. Device Parts

① PRODUCT

② USB CABLE

③ INSTRUCTIONS

4.The LED Light Therapy Mask Description

4.1 The LED Light Therapy Mask Consists Of


1. Silicone flexible face mask (Includes wireless controller + battery box)
2. USB cable and adapter
3. User manual
4. Storage bag

4.2 Features

1. Top industrial aesthetics and simple fashion design, 3D three-dimensional design fits the face better
2. Eye with medical special silicone eye pads, protect your eyes from light stimulation.
3. You can choose the appropriate working mode and adjust the time according to your skin care needs.
4. Compact and lightweight, easy to carry, so you can take care of your skin anytime, anywhere.

4.3 Remote Wireless Controller Marking

The wireless controller switches the LED ON/OFF and controls power to the mask. The wireless controller contains a rechargeable.

Sign	Meaning
ON / OFF 	<ul style="list-style-type: none"> - Long press to turn on/off, any one mode will work for 10 minutes and then turn off automatically. - Short press to switch mode, power on default 01 mode.
Light level Adjustment 	Power on the default is one, the highest is three, “-” to the first light brightness can not be decreasing, “+” to the third light brightness can not be increasing.

5. Contraindications, Warnings And Precautions

5.1 Contraindications

- Do not use the mask to treat any other conditions apart from those listed in the indications for use. The mask has not been tested for any other conditions than those listed, and the risk is unknown.
- Do not use the mask if you are pregnant, planning on becoming pregnant or breastfeeding. The mask has NOT been tested on pregnant or breast-feeding women and therefore the risk to the foetus, neonate or pregnant women is unknown.
- Do not use the mask if you suffer from Lupus erythematosus, photosensitive eczema or Albinism. If you use the mask to treat Lupus erythematosus, Photosensitive eczema or Albinism you may cause a severe skin reaction.
- Do not use the mask if you suffer from any photosensitive disorder (sensitization to light). If you use the mask and you suffer from a photosensitive disorder you may cause a severe skin reaction.
- Do not use the mask if you are taking any medication that can cause photosensitivity. If you use the mask and you are taking any medication that can cause photosensitivity you may cause a severe skin reaction. Photosensitivity is a common side effect of various medications. These can include certain antibiotics, chemotherapy drugs, and diuretics, If you are unsure about any medication you may be taking consult your healthcare provider.
- People with wounds should not use it or avoid the wound to use.
- The patient is an intended operator. The patient can safely use skin care function. The patient should not perform maintenance while the product is in use, including charging battery. Face mask is applied part.

5.2 Warnings

- Before using the mask you must carry out a sensitization test (see the following). Do not use the mask if you are sensitive to mask.
- Do not use the mask in any way that is not exactly as specified in these Instructions for use. Use of the mask other than exactly as specified in these instructions for use has not been tested and the risk is unknown.
- Turn off and remove the mask immediately if it feels hot or your skin feels uncomfortably warm.
- Do not use the mask without removing creams, lotions or make-up before treatment. Some ingredients in these products may cause light sensitization or may reduce the effect of the treatment. If you use the mask without removing these products, you may cause a skin reaction, or you may reduce the effect of the treatment.
- Never use mask near water or where water can be spilt onto the equipment
- Do not use the mask if you think it has come into contact with water.
- Keep the mask away from hot surfaces.
- Do not attempt to modify the mask in any way. Modifying the mask will void your warranty. Modifying the mask could cause an electric shock when the system is used again. The mask contains no serviceable parts.
- Do not use the mask if the machine, controller, charging plug or cable are damaged in anyway.
- Do not use the mask during charging.
- Do not touch the mask, USB cable and your USB adapter with wet hands during charging.
- Do not store flammable and explosive materials (sprayers,oil,gas,etc.)near the product.
- Do not use this product in a bath, bathroom, swimming pool, etc. where it can be directly exposed to water or liquids.
- Do not use heaters in the accessories of the product.
- Do not use the product close to or stacked with other equipment. If it must be used close or stacked, it should be observed to verify that it can operate normally under the configuration used.
- Do not place it in the reach of pets and children to avoid abusing the product or causing damage to the product.

- Do not expose the product to the outdoors for a long time to avoid contamination of the product by pests.
- Do not excessively bend or squeeze the product and power cable with heavy objects.
- Do not modify or replace the power cable and battery by yourself.
- Do not look directly at the light source part of the LED light therapy mask.
- Do not expose to harmful optical radiation.
- Do not touch cables and hoses at will to avoid strangulation caused by too long. Avoid inhalation or ingestion of small parts. Persons with potential allergies to the accessible materials used in the equipment should not use the equipment.
- Do not touch the device with wet hands when the device is in use to avoid contact injury.
- Do not use accessories removable parts and materials not described in the instructions
- Do not interconnect this device with other devices not described in the instructions for use.
- Caution-Use of controls or adjustments or performance of procedures other than those specified herein result in HAZARDOUS radiation exposure.
- Do not leave the equipment unattended when it is switched on to avoid the risk of fire or burns.
- This equipment is not intended for use by persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, instruction concerning use of the equipment to avoid the risk of fire or burns.
- This equipment is not intended for use by children. Children should be supervised to ensure that they do not play with the equipment to avoid the risk of fire and burns.
- If the mains cord is damaged, you shall have it replaced by MANUFACTURER, a service centre authorised by MANUFACTURER or similarly qualified persons in order to avoid risk of electrocution.
- Check if the voltage indicated on the equipment corresponds to the local mains voltage before you connect the equipment to avoid risk of electrocution or permanent damage to the equipment.
- If the wall socket used to power the equipment has poor connections, the plug of the equipment becomes hot. Make sure you plug the equipment into a properly installed wall socket to avoid the risk of fire and burns.
- Do not subject the equipment to heavy shocks to avoid risk of damage to the lamp.

- Always unplug the appliance after use and in case of a power failure to avoid risk of damage to the equipment.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance.
- Cleaning and user maintenance shall not be made by children without supervision.

5.3 Precautions

- Except for the connecting wires and charging wires sold by our company as spare parts of internal components, the use of accessories and cables outside the regulations may result in an increase in the emission of the LED light therapy mask.
- The purchaser or user must install and use it according to the information provided in the instruction manual.
- The product is a personal item and it is not recommended to share it with others.
- The LED light therapy mask strap contains a mixture of Polyester and nylon. If you are allergic to any of these materials, you may not be able to use the head strap provided.
- The expected operator or user must be trained in relevant skills and knowledge.
- The electric shock protection level of the remote control is Class 2, and the electric shock protection level of the mask belongs to the discharge of internal power supply of the medical device.

If device is intended for connection to a separate power supply, either the power supply shall be specified as part of the device or the combination shall be specified as a system. Users and patients/client should pay attention to eye and skin protection and do not look directly at the light when using it.

The light therapy skin care function is safe for the patient to use, and where applicable, the patient cannot safely use other unrelated functions.

Patients should not perform maintenance on their own, including changing batteries. Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

One battery cycle means that the battery is fully charged to power off.

6. Operating Mode Information

Mode of operation: Continuous operation recommended.

This product has four operating modes: Radiance Mode, Purity Mode, Deep Repair Mode, and Ultimate Revival Mode.

1. Radiance Mode with red light: wavelength 630 nm + 830 nm

Design to help stimulate the production of collagen and elastin. Intended to help fade the appearance of wrinkles and improve skin texture.

2. Purity Mode with blue light: wavelength 415 nm + 830 nm

Designed to assist in the reduction of mild to moderate inflammatory acne.

Supports the recovery of damaged tissue and helps minimize the appearance of scarring.

3. Deep Repair Mode with dark red light: wavelength 630 nm + 660 nm + 850 nm

Aims to fade the appearance of deep wrinkles and stimulate collagen and elastin production. Intended for a firming effect, reducing the look of fine lines and improving overall skin quality.

4. Ultimate Revival Mode with purple light: wavelength 415 nm + 660 nm + 830 nm

An enhanced version of Purity Mode and a combination of features from all modes. Designed to target bacteria, reduce the appearance of inflammation and redness, help prevent scarring, and support tissue regeneration.

7. Sensitization Test

Before carrying out your first treatment you must carry out a sensitization test. Follow the step by step guide below:

1. Place the machine over your inner arm.
2. Switch ON the controller and allow the machine to run for its full 10-minute treatment program.
3. At any time if your skin feels HOT or uncomfortable, STOP using the device immediately.
4. After the treatment wait at least 6 hours. If the skin on your inner arm appears red and/or itchy then LED light therapy mask is unsuitable for you. Return the device to the retailer within 14 days of purchase.


8. Using LED Light Therapy Mask

8.1 Treatment Recommendations

For best results: 10 minutes for each treatment with 3-4 treatment a week, reduce to 1-2 treatment a week once the results shown.

8.2 Carrying Out Treatment

1. Assemble your mask as described above.
2. Clean and thoroughly dry your face and hands.
3. Adjust the elastic strap on the mask to ensure the mask fits comfortably.
4. Sit or lie down and place the mask over your face.
5. Use the controller to select light color, brightness, flickering or steady light and treatment time.
6. Apply your preferred skincare creams after treatment as absorption is enhanced following light therapy.
7. Continue with your skin care routine. 10 minutes for each treatment with 3-4 treatment a week, reduce to 1-2 treatment a week once the results shown.

8.3 Tips

If you use the mask in manual mode, it should be controlled within 10 minutes. If it is used for too long, it may cause damage to the skin.

If you are using the product for the first time, please refer to the following maintenance itinerary and adjust it according to your skin condition.

10 minutes for each treatment with 3-4 treatment a week, reduce to 1-2 treatment a week once the results shown.

9. Using The Wireless Controller

9.1 Modes

The photomask is divided into 4 modes, long press ON/OFF, the power on defaults to Radiance mode, short press ON/OFF to switch the mode. Any of the modes will be turned off automatically after 10 minutes of operation.

9.2 Light Level Selection

There are 3 levels to adjust the light intensity, the default is 1 level, “ - ” to 1 level cannot be reduced, “ + ” to 3 levels can not increase the brightness of the light.

9.3 Voice Announcement

The default language of the mask is French. To switch to English, press and hold the “-” button for 5 seconds. To return to French, repeat the same action.

-Voice message when powered on: Welcome to your LED Face Mask. Radiance mode activated.

-Mode 02, voice message: Purity Mode activated

-Mode 03, voice message: Deep Repair Mode activated

-Mode 04, voice message: Ultimate Revival Mode activated

-Mode 01, voice message: Radiance Mode activated

-Voice message for the first light level: Level 1

-Second level: Level 2


-Third level: Level 3

-When the battery is low: Low battery, please charge.

-When powering down : Powering down.

10. Charging

Warm Reminder: Before charging, prepare the mask, USB cable.

There is a battery box on the mask, connect the magnetic tip of the equipped USB cable to the charging port of the battery box, connect the port of the USB to the adapter, you need to make sure that the voltage labeled on the product is the same as the power supply voltage. The mask cannot be activated during charging.

Breathing light:

- Under power-on state, the breathing light presents the state as a sense of breathing.

- Under charging state, the breathing light presents a sense of breathing.

- When the battery is fully charged, the breathing light will be on for a long time. Low voltage alarm:

- When the battery voltage is lower than 3.3V during power-on use, the voice broadcasts three times to indicate low voltage, 3.1V voice broadcasts shutdown, automatic shutdown, and the breathing light flashes five times.

- After the low voltage shutdown, if the charging operation is not carried out it is impossible to start the phototherapy cover again.

Battery type: lithium ion battery.

WARNING : Avoid using this device next to or stacked on top of other devices as this may result in improper operation.

If this equipment must be used, observe this and other equipment to verify that they are operating properly”

WARNING : Use of accessories, sensors, and cables other than those specified or supplied by the manufacturer of this equipment may result in misoperation.

WARNING : Use of accessories, sensors, and cables other than those specified or supplied by the manufacturer of this equipment may result in malfunction.

WARNING : Use of accessories, sensors, and cables other than those specified or supplied by the manufacturer of this equipment may result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and may result in mishandling.”

WARNING : When using portable RF communication equipment (including peripheral equipment such as antenna cables and external antennas), it must not be located closer than 30 cm (12 inches) from any part of the ME device, including cables specified by the manufacturer.

Failure to do so may result in degraded performance of this equipment.

Table 1

Declaration-electromagnetic emission	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emission- sIEC61000-3-2:2009	Class A
Emissions test	Clause 5

Table 2

Declaration-electromagnetic emission		
Immunity test	IEC 60601 test level	Compliance level

Electrostatic discharge (ESD) IEC61000-4-2:2008	±8 kV contact ±2 kV,±4 kV,±8 kV, ±15 kV air	±8 kV contact ±2 kV,±4 kV,±8 kV ±15 kV air
Electrical fast transient/burst IEC61000-4-4:2012	±2 kV for power supply lines ±1 kV for input/out- put lines	±2 kV for power supply lines ±1 kV for input/out- put lines
Surge IEC61000-4-5:2014	±0.5kV, ±1kV line(s)to lines±0.5kV, ±1kV, ±2 kV line(s)to earth	±0.5kV, ±1kV line(s)to lines±0.5kV, ±1kV ±2 kV line(s)to earth
Voltage dips, short interruptions and Vo–ltage variations on power supply input lines IEC61000-4-11:2004	0%UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0%UT;1cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 250/300 cycles	±8 kV contact ±2 kV,±4kV,±8 kV, ±15 kV air
Power frequency (50/60 Hz)magnetic field IEC61000-4-8:2001	30 A/m	30 A/m
NOTE:UT is the a.c.mainsvoltage prior to application of the test level.		

Table 3

Declaration–electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level

IConducted RF IEC61000-4-6:2006	3V 3V 0.15 MHz to 80 MHz 6 V in ISM bands between 0.15 MHz and 80MHz	3V 0.15 MHz to 80 MHz 6V in ISM bands between 0.15 MHz and 80 MHz
Radiated RFIEC61000-4-3:2010	10V/m 80 MHz to 2.7GHz	10 V/m

Table 4

Declaration-IMMUNITY to proximity fields from RF wireless communications equipment					
Immunity test	IEC60601 test level				Compliance level
	Test frequency	Modulation	Maxim power	Immunity level	
Radiated RFIEC 61000- 4-3:2010	385 MHz	**Pulse Modulation: 18Hz	1,8 W	27 V/m	27 V/m
	450 MHz	*FM+5Hz deviation: 1kHz sine	2 L	28 V/m	28 V/m
	710 MHz 745 MHz 780 MHz	**Pulse Modulation: 217Hz	0,2 W	9 V/m	9 V/m
	810 MHz 870 MHz 930 MHz	**Pulse Modulation: 18Hz	2 L	28 V/m	28 V/m
		**Pulse Modulation: 217Hz	2 L	28 V/m	28 V/m

Note*-As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case. Note**-The carrier shall be modulated using a 50% duty cycle square wave signal.

Tips:

- When charging, you only can use the original USB cable equipped with the mask and the mask should not be used during charging.
- Since the battery capacity of this product is 1800mAh, we recommend you to use the adapter included or a 0.5A adapter.
0.5A adapter to reduce charging time
- When the battery voltage is lower than 3.2V, the indicator light on the remote control will flash and the mask will automatically shutdown after flashing 10 times.

11. Risks And Benefits

Known Risks

- Immediately after treatment skin may appear red. Skin should return to normal within 1-2 hours.

Known Benefits

- Reduce facial fine lines and improve texture and pigmentation.
- Improve skin imperfection.

12. Safety Signs And Meaning

Safety sign	Meaning
	Electrical appliances for indoor use.
	Don't throw away at will.
	Compliant with CE's EMC standard requirements
	This symbol on the device label means: Read the user
	This symbol means keep away fire.
	This symbol means place upwards.

	<p>This symbol means keep dry.</p>
	<p>TYPE BF APPLIED PARTS</p>
<p>IP22</p>	<p>IP classification</p>

13. Care And Maintenance

- Clean the surface of the LED light therapy mask after each treatment
- Store the LED light therapy mask at room temperature in dry conditions
- Avoid storing in dusty conditions, indirect sunlight or in damp and humid conditions.
- The LED light therapy mask contains no serviceable parts.
- Before each use, with the device switched off, clean the face mask following the cleaning instructions below.

13.1 Cleaning Of LED Light Therapy Mask

- (1) Please turn off the power when cleaning and tidying up the product.
- (2) Use a 70%-80% ethanol solution (rubbing alcohol) to wipe the surface of the mask twice for 3 minutes, and then wipe with a damp or dry cloth.
- (3) When cleaning this product, do not use benzene, thinner, butyl acetate and other flammable materials to wipe.
- (4) When cleaning the product, do not spray water directly on the product. First, use a 70%-80% ethanol solution to wipe the silica gel surface twice for 3 minutes, and then wipe it with a damp or dry cloth.
- (5) DO not use any other cleansing products as these may damage the plastics.

13.2 Disposal

After the service life of the equipment expires, if the materials composed of the lithium battery and the instrument are discarded everywhere, it will pollute the environment. Regarding the disposal of packaging materials and end-of-life products, remember to follow local laws and regulations. Do not dispose of any

part of the machine at a domestic waste facility. The laws and ordinances based on users require proper treatment of end-of-life products and materials and try to support their sorting and recycling work.

14. Troubleshooting

Problem	Possible cause	Action
Can't turn on the power	<ol style="list-style-type: none"> 1. The power plug is not properly inserted. 2. The LED light therapy mask is not connected properly. 3. The power button was not long pressed for one second. 	<ol style="list-style-type: none"> 1. Please insert the power plug correctly. 2. Please check the power cable. 3. Please long press the power button again. 4. Please charge the remote control.
LED light does not light up	<ol style="list-style-type: none"> 1. The LED light therapy mask is not connected properly 2. The circuit board is malfunctioning. 3. Pull the power cord excessively 	<ol style="list-style-type: none"> 1. Please check the power cable 2. Do not pull the power cable excessively.
Flashing appears during use. Cannot switch the lighting mode.	<ol style="list-style-type: none"> 1. The battery voltage is too low of the remote control. 2. The power cable is not fully inserted into the remote controlport. 3. Device is faulty. 	<ol style="list-style-type: none"> 1. Please charge the remote control. 2. Please check if the power cable is in good condition. 3. Return to your supplier.

If the problem reoccurs, and your device is within the warranty period please return to your provider.

- DO NOT attempt to investigate any further.
- DO NOT attempt to modify the device in any way. The LED System has NO parts you can repair yourself.

15. Environmental Conditions

Item		
Operation environment	Temperature	5°C~+35°C
	Humidity	85%RH or less
	Atmospheric pressure	80-106 KPA
	Altitude	≤2000 meters
Storage environment	Temperature	-10°C~+50°C
	Humidity	85%RH or less
	Atmospheric pressure	80-106 KPA
	Altitude	≤2000 meters

The external environment that the product can withstand

(1) The mask can be placed in an environment with a temperature of $60\pm 2^{\circ}\text{C}$ and a humidity of $85\%\pm 2\%\text{RH}$ for 48 hours, and then placed at room temperature for 2 hours.

(2) The mask can be placed in an environment with a temperature of $-20\pm 2^{\circ}\text{C}$ for 48 hours, and then placed at room temperature for 2 hours.

(3) The mask can be placed in an environment with a temperature of $40\pm 3^{\circ}\text{C}$ and a humidity of $90\%\pm 2\%\text{RH}$ for 48 hours, it is placed at room temperature for 2 hours.

(4) The mask can be placed in an environment with a temperature of 40°C for 2 hours, and test related functions in this environment.

(5) The mask can be placed in an environment with a temperature of 0°C for 2 hours, and test related functions in this environment.

16. EMC Manufacturer's Declaration And Guidance

1. The mask needs to take special precautions for electromagnetic compatibility (EMC), and the purchaser or user must install and use it according to the information about EMC compliance (YY0505-2012) provided in the instruction manual.

2. Portable and mobile radio frequency communication equipment may affect the facial photon therapy device. Please keep away from such equipment when using it.

3. In addition to the connecting wires and charging cables sold by our company as spare parts of internal components, the use of accessories and cables outside the regulations may result in an increase in the emission of the facial photon therapy device or a decrease in immunity.

The following cable types must be used to ensure compliance with and adherence to interference radiation and immunity standards:

Cable	Length(m)
Power cord	0,5
Charging Cable	0,5

4. This mask should not be used close to or stacked with other equipment. If it must be used close or stacked, it should be observed to verify that it can operate normally under the configuration used.

5. This mask should not touch the pins of the connectors marked with the electrostatic discharge warning symbol, and should not be connected to these connectors unless electrostatic discharge precautions are used.

Provisions on preventive measures against electrostatic discharge:

Human body or objects carry different electrostatic voltages due to charge transfer; since electrostatic discharge is completed in the order of ns or us, the peak current can reach tens of amperes, and the instantaneous power is very huge. The resulting electrostatic discharge electromagnetic pulse energy is enough to damage the sensitive components in the electronic components; because the rise time of the current waveform is very short, that is, the rate of change of the current (di/dt) is very large, it can induce a high potential of hundreds of volts or even thousands of volts, thereby a strong electric field is generated to breakdown the sensitive element. In order to prevent damage to the instrument, the following measures should be taken:

- (1) Ensure the environmental humidity;
- (2) Laying anti-static floor or carpet;
- (3) The operator should wear an anti-static hand strap on the wrist. This hand strap should have good grounding performance. This measure is the most effective.

6. The basic performance of the facial photon therapy instrument: red light mode: radiation intensity $E_e \leq 96\text{mW/cm}^2$. blue light mode: radiation intensity: $E_e \leq 120\text{mW/cm}^2$.

7. The guidelines and the company's statement are shown in Table 1 to Table 4; This mask meets the electricity magnetic compatibility requirements of IEC 60601-1-2 /AMD 1:2020 PRV idtIEC 60601-1-2:2004; Guidelines and manufacturer's declaration-electromagnetic emissions-for all equipment and systems.

Table 1 Electromagnetic radiation
(corresponding to Table 201 in IEC 60601-1-2/AMD1:2020 PRV)

Guidelines and manufacturer's declaration-electromagnetic emissions		
The mask is expected to be used in the electromagnetic environment specified below, The purchaser or user of the face photon therapy device should ensure that it is used in the electromagnetic environment.		
Launch test	Compliance	Electromagnetic environment-guide
RF emission CISRP11:2016	Team 1	The facial photon therapy device uses radio frequency energy only for its internal functions. Therefore, its radio frequency emission very low, and the possibility of causing interference to nearby electronic equipment is very small.
RF émission CISRP11:2016	Type B	The facial photon therapy device is suitable for use in all facilities, including households and directly connected to the residential
Harmonic emission IEC61000-3-2:2009	Type A	
Voltage fluctuation/ flicker emission IEC61000-3-2:2009	Correspond	Public low-voltage power supply network for households.

Guidelines and manufacturer's declaration–electromagnetic immunity–for all equipment and systems, this mask is expected to be used in the electromagnetic environment specified below. The purchaser or user of this mask should ensure that it is used in this electromagnetic environment.

Table 2 Electromagnetic immunity 1
(corresponding to table 202 in IEC60601-1-2/AMD1:2020 PRV)

Guidelines and manufacturer's declaration–electromagnetic immunity			
The mask is expected to be used in the electromagnetic environment specified below, and the purchaser or user of the face photon therapy device should ensure that it is used in this environment.			
Immunity test	IEC60601 test level	Coincidence level	Electromagnetic environment-guide
Electrostatic discharge IEC 61000-4-2:2008	±6kV Contact discharge ±8kV Air discharge	±6kV Contact discharge ±8kV Air discharge	The floor should be wood, concrete, or ceramic tiles. If the floor is covered with a contaminated material, the relative humidity should be at least 30%.
Voltage transient burst IEC 61000-4-4:2012	±2kV To the power cord ±1kV To input/output cord	±2kV To the power cord Not applicable	The network power supply should have the quality used in a typical commercial or hospital environment.
Surge IEC 61000.4-5:2014	±2kV line to line ±1kV line to ground	±2kV line to line Not applicable	The network power supply should have the quality used in a typical commercial or hospital environment.

<p>The voltage on the voltage input line is temporarily low, short-term interruption and voltage change IEC 61000-4-11:2004</p>	<p><5%UT, lasting for 0.5 weeks (at UT>95% temporary low) 40% UT for 5 weeks(on UT, a temporary low of>60%)70% UT, for 25 weeks (at UT, a temporary low of>30%) 5%UT, lasting 5s (at UT atemporary</p>	<p><5%UT, lasting for 0.5 weeks (at UT>95% temporary low) 40% UT for 5 weeks(on UT, a temporary low of>60%)70% UT, for 25 weeks (at UT, a temporary low of>30%) 5%UT, lasting 5s (at UT atemporary</p>	<p>The network power supply should have the quality used in a typical commercial or hospital environment. If the user of the facial photon therapy device needs to continue to operate during the power interruption, it is recommended that the facial photon therapy device be powered by an uninterrupted power supply or battery.</p>
<p>Power frequency magnetic field (50Hz) IEC 61000-4-8:2001</p>	<p>3A/m</p>	<p>3A/m</p>	<p>If image distortion occurs, the facial photon therapy device must be placed farther away from the source of the power frequency magnetic field or a magnetic shielding protection must be installed. The power frequency magnetic field should be measured to ensure that the installation position is sufficiently low when it is expected to be used</p>
<p>Note: UT refers to the AC network voltage before the test voltage is applied</p>			

Guidance and manufacturing statement-electromagnetic immunity-for non-life support equipment and systems.

Table 3 Electromagnetic immunity (corresponding to Table 204 in IEC 60601-1-2/AMD1:2020 PRV)

<p>Guidelines and manufacturer's declaration-electromagnetic immunity</p>
<p>The mask is expected to be used in the electromagnetic environment specified below, and the purchaser or user of the mask should ensure that it is used in this environment.</p>

Immunity test	IEC60601 test level	Coincidence level	Electromagnetic environment-guide
Radio frequency conduction IEC 61000-4-6:2006	3V(Effective value) 150KHz-80MHz	3V (Effective value)	<p>Portable and mobile radio frequency communication equipment should not be used closer to any part of the facial photon therapy device, including cables, than the recommended isolation distance. This distance should be calculated based on the formula corresponding to the transmitter frequency. Recommended isolation distance $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ 80MHz-800MHz $d=2.3\sqrt{P}$ 80MHz-2.5GHz In the formula: P-According to the maximum rated output power of the transmitter provided by the transmitter manufacturer, in watts; d-Recommended isolation distance, in meters, the field strength of the fixed radio frequency transmitter is determined by surveying</p>
<p>Note 1: At 80MHz and 800MHz frequency points, the higher frequency band formula is used Note 2: These guidelines may not be suitable for all situations. Electromagnetic propagation is affected by absorption and emission from buildings, objects and humans.</p>			
<p>a:Fixed transmitters, such as base stations for wireless (cellular/cordless) telephones and terrestrial mobile radios, business radios, AM and FM radio broadcasting, and TV broadcasting. The field strength cannot be accurately predicted theoretically. In order to assess the electromagnetic environment of fixed radio frequency transmitters, the survey of electromagnetic fields should be considered. If the measured field strength of the place where the facial photon rejuvenation device is located is higher than the above-mentioned applicable radio frequency compliance level, the facial photon treatment device should be observed to verify its normal operation. If abnormal performance is observed, supplementary measures may be necessary, such as re-adjusting the direction or position of the facial photon therapy device. b: In the entire frequency range of 150kHz-80MHz, the field strength should be less than 3V/m</p>			

Recommended isolation distance between portable and mobile RF communication equipment and this mask or system-for non-life support equipment or systems.

Table 4 Safety distance
(corresponding to table 206in IEC 60601-1-2/AMD1:2020 PRV)

Recommended isolation distance between portable and mobile RF communication equipment and facial photon therapy device or system			
The mask is expected to be used in an electromagnetic environment where radiated RF disturbances are controlled. According to the maximum output function of the communication device, the purchaser or user of the facial photon therapy device can prevent electromagnetic interference by maintaining the minimum distance			
Electrostatic discharge IEC 610004-2:2008	Corresponding to the isolation distance of the transmitter at different frequencies/m		
	150 kHz à 80 MHz d=1,2 Vp	80 MHz à 800 MHz d=1,2 Vp	80 MHz à 800 MHz d=2,3 Vp
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
0.01	1.2	1.2	2.3
0.01	3.8	3.8	7.3
100	12	12	23

For the maximum rated output rate of the transmitter not listed in the above table, the recommended isolation distance d , in meters(m), can be determined by the formula in the corresponding transmitter frequency column, where P is provided by the transmitter manufacturer. The maximum rated output power of the transmitter, in watts (W).

Note 1: At 80MHz and 800MHz frequency points, the formula for the higher frequency range should be used.

Note 2: These guidelines may not be suitable for all situations. Electromagnetic propagation is affected by absorption and reflection from buildings, objects, and humans.

17. Service Life

The service life of the product is three (3) years. See FAQ for more details.

18. Guarantee Information

1. Our company guarantees that the indicators of the products will be strictly inspected when the products leave the factory to ensure that the qualification rate of the products reaches 100%.
2. Our company guarantees that within one year from the date of purchase, if the product itself has quality problems, our company will repair it free of charge.
3. Our company can provide users with the necessary information so that Our company can provide users with the necessary information so that qualified technicians can carry out maintenance on the parts of the equipment designated by our company as repairable.
4. The charging cable can be replaced by yourself, users can buy the quality-assured 5V1A magnetic charging cable from regular manufacturers.
5. If there is a problem with the wireless remote control, do not replace it yourself. If replacement is necessary, please contact the manufacturer.
6. Warranty Exclusion (International Use): Please note that this warranty does not cover any damage to the battery or electronic components resulting from charging the device outside of Canada. The use of voltage converters or international plug adapters is at the user's own risk and may void the warranty.

* The warranty does not apply in the event of wear and tear, incorrect voltage usage, visible damage, or other misuse.

19. Frequently asked questions

1. What does the recommended usage look like?

For best results, use the mask for a 10-minute session, 3 to 5 times per week. After that, reduce to 1 or 2 sessions per week, or as needed, to maintain your results.

2. How soon can results be seen?

The appearance of results varies from person to person – generally between 10 and 30 days, depending on consistency and frequency of use. Optimal results typically appear after several weeks, often around 8 weeks of consistent use. To maintain the benefits, it's important to continue using the mask regularly.

3. What are the wavelengths for each mode?

- Radiance Mode (Red Light): 630 nm + 830 nm
- Purity Mode (Blue Light): 415 nm + 830 nm
- Deep Repair Mode (Deep Red Light): 630 nm + 660 nm + 850 nm
- Ultimate Revival Mode (Purple Light): 415 nm + 660 nm + 830 nm

4. What does “nm” mean?

“Nm” stands for nanometer, a unit used to measure the wavelength of light. Depending on the wavelength, the light penetrates to different depths of the skin, acting on various layers and addressing specific skin concerns.

5. How many LEDs does the mask contain?

The mask contains 60 LED units, for a total of 240 light bulbs.

6. What is the irradiance in mW/cm²?

It varies depending on the selected mode:

- 630 nm + 830 nm = 31.4 mW/cm² ± 25%
- 415 nm + 830 nm = 43.4 mW/cm² ± 25%
- 630 nm + 660 nm + 830 nm = 57.2 mW/cm² ± 25%
- 415 nm + 660 nm + 830 nm = 78.6 mW/cm² ± 25%

This power level allows for excellent light penetration while remaining safe and effective.

7. Can I use it if I am photosensitive?

No, to avoid any skin reactions.

8. Sensitive skin: is there a risk of reaction?

It is recommended to perform a sensitivity test on the inside of your arm before the first treatment. For your first session, start at the lowest intensity and observe your skin's reaction. If no irritation occurs, you can gradually increase the frequency and/or intensity.

9. Is it recommended for someone who is 20 years old?

The minimum recommended age is approximately 18 years, as the skin is mature enough for this type of treatment. For teenagers struggling with acne, the Purity Mode may be used under supervision.

10. I am pregnant, can I use it? What about breastfeeding?

It is always best to consult a healthcare professional before using the mask. Red light therapy is generally considered safe during pregnancy if there are no specific skin conditions.

11. Is the mask suitable for rosacea? Can it replace IPL treatment?

LED technology, specifically red and infrared light, helps soothe the skin and reduce the appearance of redness and discomfort. While not a medical treatment for rosacea or a substitute for clinical IPL, this mask is an ideal at-home aesthetic tool to promote a more uniform, decongested-looking complexion.

12. What is the battery life?

The battery life is approximately 4–6 uses per charge at light intensity level 3, and 15 uses at light intensity level 1. The battery lasts for over 500 full charging cycles. We recommend waiting until the battery is low before recharging the mask.

13. What does the service life mean?

The three-year useful lifespan corresponds to the period during which the mask delivers optimal performance in terms of light intensity, LED efficiency, and battery capacity, based on your usage. This does not mean that it stops working after three years, but rather that this is the timeframe during which its maximum effectiveness is guaranteed. The actual lifespan may vary depending on usage frequency and maintenance conditions.

14. Can I travel with my LED mask?

Yes, the LED mask is designed for international travel. Although the device supports a voltage range of 100V to 240V, it is important to use an appropriate voltage converter and local plug adapter to ensure optimal performance and protect the battery, as you would with any other electrical device.

15. Does the LED mask help with melasma?

LED light may help improve the appearance of melasma by promoting cellular regeneration and evening out skin tone. However, it does not replace a targeted medical treatment. Since melasma is a complex and often hormonal condition, it is always best to consult a dermatologist or healthcare professional to determine the most appropriate treatment plan.

Made in China

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1650 avenue Jules-Vernes, Local 103. G2G 2R1, Québec, Canada